



HPLC and UHPLC Methods for Melamine



Milk, infant formula and other dairy products were recently found contaminated with melamine, following an earlier melamine contamination outbreak in pet food in 2007.


Toxicology studies show ingestion of melamine in large quantities may lead to reproductive damage or bladder cancer due to the formation of bladder/kidney stones. At lower levels, melamine and cyanuric acid are absorbed into the bloodstream. Together, they concentrate and form melamine cyanurate in the urine-filled renal microtubules. Crystallization blocks and damages the renal cells that line the tubes, causing the kidneys to malfunction.

After 2007 and more recent melamine contamination outbreaks, there is an urgent need for analytical methods that can identify and quantify melamine in food. Current melamine detection methods involve LC-MS and GC-MS. GC-MS methods require derivatization, and LC-MS methods generally use gradient conditions that require column clean up and re-equilibration.


Grace has developed MS-compatible LC methods for Melamine using an HILIC media platform that can be applied to both traditional HPLC as well as UHPLC systems. Melamine was analyzed with a standard HPLC system using a 5µm particle HILIC phase packed into a 250 x 4.6mm column. The 1.5µm version of this phase was then packed into a high throughput format conducive to UHPLC and fast LC systems. Both methods deliver excellent linearity and use isocratic elution for fast analysis without the need for re-equilibration.

HPLC Method for Melamine

This HPLC analytical method for melamine fulfills the FDA requirements using a HILIC column and an ionizable mobile phase compatible with mass spec. Low UV detections offers excellent linearity between 40ng and 100µg.




HPLC Column: VisionHT™ HILIC, 5µm, 4.6 x 250mm (Part No. 86466)
Mobile Phase: Acetonitrile:10mM Ammonium Acetate in Water (95:5)
Flow Rate: 1mL/min
Detection: UV@240nm
Column Temperature: 30°C
Injection: 40µg/mL x 20µL




UHPLC Method for Melamine

Compared to the conventional HPLC method, the UHPLC method is 4 times faster. With the use of 1.5µm particles, optimal linear velocities extend over a wider range. Therefore, it is possible to maintain efficiency and resolution while running samples at faster flow rates.




UHPLC Column: VisionHT™ HILIC, 1.5µm, 2 x 50mm (Part No. 5141919)
Mobile Phase: Acetonitrile:10mM Ammonia Acetate in Water (95:5)
Flow Rate: 0.2mL/min
Detection: UV@240nm
Column Temperature: 30°C
Injection: 50µg/mL x 0.5µL

9 injections in parallel shows good reproducibility



UHPLC Column: VisionHT™ HILIC, 1.5µm, 2 x 50mm (Part No. 5141919)
Mobile Phase: Acetonitrile:Water(20mM Ammonium Formate) (90:10)
Flow Rate: 0.2mL/min
Detection: UV@240nm
Column Temperature: 30°C
Injection: 50µg/mL x 0.1µL




Conc.	Inj. (µL)	Weight (ng)	Peak Area
50µg/mL	0.1	5	7.3
50µg/mL	0.5	25	32
50µg/mL	1	50	63
50µg/mL	5	250	318
50µg/mL	10	500	641
50µg/mL	20	1000	1293

This method exhibits excellent linear response between 5ng and 1000ng for accurate quantitation.


HPLC Method Transfer to UHPLC

The recent adoption of Fast LC systems within the laboratory environment means there is a period of time when most labs have a variety of LC systems types – Traditional, UHPLC, and/or alternate Fast LC systems. Optimizing and transferring methods between systems has not been a simple and intuitive task. However, when the identical bonded phase is available in sub2, 3, 5, and 10µm particle sizes, it can be applied to appropriate formats to suit the system type. A simple calculation to determine equivalent linear velocity is all that's necessary to seamlessly transfer methods between systems, or to optimize a method on one system and apply to an entirely different type of system. Compared to a conventional HPLC method, the UHPLC method is 4 times faster.



	HILIC Column	Flow Rate	Time (min)	Conc.	Inj. Range	Loading Range
HPLC	VisionHT™ HILIC, 5µm, 4.6 x 250mmL	1.0mL/min	12.090	40µg/mL	1-100µL	40ng-4mg
UHPLC	VisionHT™ HILIC, 1.5µm, 2.0 x 50mmL	0.2mL/min	3.259	50µg/mL	0.1-20µL	5ng-100ng

ELSD Methods for Melamine



HPLC Column: VisionHT™ HILIC, 5µm, 4.6 x 250mm (Part No. 86466)
Detector: Alltech® 3300 ELSD (Part No. 5135834)
Mobile Phase: Acetonitrile:10mM Ammonium Acetate in Water (95:5)
Flow Rate: 1mL/min
Detection: UV@240nm
Column Temperature: 30°C
ELSD 3300 Settings: drift tube 40°C, gas 1.8L/min, gain x 4
Injection: 40µg/mL x 20µL

SPE Methods for Melamine

Extracting melamine from various food matrices using solid phase extraction is often required to obtain accurate quantitation. GracePure™ Cation-X, a strong cation exchange resin delivers a cleaner, more concentrated sample.


Melamine from Ham

SPE Column Type: GracePure™ SPE Cation-X, 500mg/3mL (Part No. 5138770)

Sample Pretreatment: weight ham 2g
 mix with 20mL 1% Trichloroacetic Acid and 5mL Acetonitrile
 Sonicate for 30mins
 Centrifuge for 15mins
 Filter the supernatant with 0.45µm syringe filter, then spike with melamine standard solution

SPE Steps:

Column Conditioning: 2 column volumes methanol, then 2 column volumes water
Sample Application: slowly apply 10mL ham extract through the column
Column Washing: 3 column volumes water, then 2 column volumes methanol
Eluent: 4mL Methanol:0.5M Ammonium Hydroxide (1V:1V)



Mobile Phase: 10mM Ammonium Acetate:Acetonitrile (5:95)
HPLC Column: VisionHT™ HILIC, 5µm, 4.6 x 250mm
Detector: UV @ 240nm
Flow Rate: 1mL/min
Column Temperature: 30°C
Injection: 10µL

SPE Methods for Melamine continued

Melamine from Biscuit

SPE Column Type: GracePure™ SPE Cation-X, 500mg/3mL, (Part No. 5138770)

Sample Pretreatment: weight biscuit 2g
mix with 20mL 1% Trichloroacetic Acid and 5mL Acetonitrile
Sonicate for 30mins
Centrifuge for 15mins
Filter the supernatant with 0.45µm syringe filter, then spike with melamine standard solution


SPE Steps:

Column Conditioning: 2 column volumes methanol, then 2 column volumes water

Sample Application: slowly apply 10mL biscuit extract through the column

Column Washing: 3 column volumes water, then 2 column volumes methanol

Eluent: 4mL Methanol:0.5M Ammonium Hydroxide (1V:1V)



SPE recovery*=80%

Mobile Phase: 10mM Ammonium Acetate:Acetonitrile (5:95)

Column: VisionHT™ HILIC, 5µm, 4.6 x 250mm

Detector: UV @ 240nm

Flow Rate: 1mL/min

Column Temperature: 30°C

Injection: 10µL

Melamine from Pork

SPE Column Type: GracePure™ SPE Cation-X, 500mg/3mL (Part No. 5138770)

Sample Pretreatment: weight pork 2g
mix with 20mL 1% Trichloroacetic Acid and 5mL Acetonitrile
Sonicate for 30mins
Centrifuge for 15mins
Filter the supernatant with 0.45µm syringe filter, then spike with melamine standard solution


SPE Steps:

Column Conditioning: 2 column volumes methanol, then 2 column volumes water

Sample Application: slowly apply 10mL pork extract through the column

Column Washing: 3 column volumes water, then 2 column volumes methanol

Eluent: 4mL Methanol:0.5M Ammonium Hydroxide (1V:1V)



SPE recovery*=103%

HPLC Mobile Phase: 10mM Ammonium Acetate:Acetonitrile (5:95)

Column: VisionHT™ HILIC, 5µm, 4.6 x 250mm

Detector: UV @ 240nm

Flow Rate: 1mL/min

Column Temperature: 30°C

Injection: 10µL

Conclusions

Fast, isocratic UHPLC and HPLC methods were developed for the determination of melamine. Linearity range for melamine is 40 - 4000 ng. Since the Grace® VisionHT™ HILIC UHPLC and the HILIC Analytical HPLC columns are based on the same media, methods developed on the HPLC column can be easily transferred to UHPLC and vice versa. Compared to the HPLC method, UHPLC is 4 times faster. Mobile phases used in these methods are MS and ELSD compatible.

www.discoverysciences.com

The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn. and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification.

ALLTECH® is a trademark, registered in the United States and/or other countries, of Alltech Associates, Inc. VISIONHT™ is a trademark of Alltech Associates, Inc. GRACE® and GRACE DAVISON® are trademarks, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. GRACEPURE™ and GRACE DAVISON DISCOVERY SCIENCES™ are trademarks of W. R. Grace & Co.-Conn. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership. Grace Davison Discovery Sciences is a product group of W. R. Grace & Co.-Conn. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. © Copyright 2009 Alltech Associates, Inc. All rights reserved. The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. Grace reserves the right to change prices and/or specifications without prior notification. Printed in the USA.

12/2009, M287

W. R. Grace & Co.-Conn. 7500 Grace Drive, Columbia, MD 21044

Grace Davison Discovery Sciences Regional Headquarters:

IN THE AMERICAS:
2051 Waukegan Road
Deerfield, IL 60015-1899
Tel: 1 847 948 8600
Fax: 1 847 948 1078
Email: discoverysciences@grace.com

IN EUROPE:
Brandstraat 12
B-9160 Lokeren, Belgium
Tel: +32 (0)9-340-65 65
Fax: +32 (0)9 340 65 60
Email: discoverysciences.EU@grace.com

IN ASIA:
19th Floor, K.Wah Center
1010 Huai Hai Zhong Road
Shanghai, 200031 PRC
Tel: 86 21 5467 4678
Fax: 86 21 5405 1500
Email: dsbiz.asia@grace.com

IN INDIA:
17 Commerce Center
Opposite Krishna Hospital, Paud Road
Pune, 411038 India
Tel: +91 20 6644 9900
Fax: +91 20 2544 1740
Email: pune@grace.com

IN AUSTRALIA/NEW ZEALAND:
2 Kerr Court
Rowville, 3178
Victoria, Australia
Tel: +61 3 9237 6100
Fax: +61 3 9237 6101
Email: discoverysciences.AU@grace.com