

GraceResolv™ High Resolution Flash Cartridges

- Higher Resolution
- Increased Loading Capacity
- Reproducible Performance
- Universal Luer Fittings




Superior Silica



Grace's 30 years of chromatography experience combined with 80 years of silica manufacturing expertise has yielded the ideal flash cartridge. GraceResolv™ cartridges are packed with a new grade of Davisil® silica designed to maximize performance of today's flash instruments.

Increase Resolution and Product Purity with GraceResolv™ Flash Cartridges



GraceResolv™ Silica Cartridges – Luer Endfittings



GraceResolv™ Cartridges – Luer Endfittings



Luer Lock Solid Loader



Luer Slip Solid Loader

GraceResolv™ Silica Cartridges – Luer Endfittings Adaptors to Connect to Biotage® Systems



Support Rings



Adaptors



Adjustable Plunger

Support Rings	
Description	Part No.
Pack contains adaptors for five sizes (4g, 12g, 40g, 80g, 120g)	8618945

Adaptors	
Description	Part No.
Male and female luer solvent line connectors	8618940

Adjustable Plunger	
Description	Part No.
3g (12mL cartridge)	8619011
25g (60mL cartridge)	8619010

GraceResolv™ Cartridges – Luer Endfittings		
Cartridge Size	Qty.	Part No.
GraceResolv™ 4g Silica Cartridge	20/pk	5154063
GraceResolv™ 12g Silica Cartridge	20/pk	5154064
GraceResolv™ 24g Silica Cartridge	15/pk	5155963
GraceResolv™ 40g Silica Cartridge	15/pk	5154065
GraceResolv™ 80g Silica Cartridge	12/pk	5154066
GraceResolv™ 120g Silica Cartridge	10/pk	5154067
GraceResolv™ 220g Silica Cartridge	6/pk	5155964
GraceResolv™ 330g Silica Cartridge	4/pk	5154068
GraceResolv™ 750g Silica Cartridge	4/pk	5155971
GraceResolv™ 1500g Silica Cartridge	3/pk	5155972

Solid Loaders				
Description	Empty*		Silica Filled	
	Qty.	Part No.	Qty.	Part No.
<i>Luer Slip Solid Loader – Isco® Compatible</i>				
3g (12mL cartridge)	100/pk	5142050	20/pk	5152371
25g (60mL cartridge)	100/pk	5142051	16/pk	5142034
<i>Luer Lock Solid Loader – Biotage® Compatible</i>				
5g (12mL cartridge)	100/pk	8618941	20/pk	8618942
25g (60mL cartridge)	100/pk	8618943	16/pk	8618944

*Empty solid loaders come with two 20µm frits

tech tip

To connect a GraceResolv™ cartridge to a Biotage® instrument, simply connect the adaptors to the inlet and outlet side of the cartridge and attach to the instrument's tubing. The cartridge then fits on the Biotage® instrument with the appropriate support ring.

If needed, a luer lock solid loader can also be used with the adjustable plunger.



www.discoverysciences.com

DAVISIL®, GRACE®, and GRACE DAVISON® are trademarks, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. GRACE DAVISON DISCOVERY SCIENCES™ and GRACERESOLV™ are trademarks of W. R. Grace & Co.-Conn. ISCO®, COMBIFLASH®, and COMPANION® are trademarks, registered in the United States and/or other countries, of Teledyne Isco, Inc. Corporation. BIOTAGE® is a trademark, registered in the United States and/or other countries, of Biotage AB. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. Grace Davison Discovery Sciences is a product group of W. R. Grace & Co.-Conn., which now includes all product lines formerly sold under the Alltech brand. © Copyright 2010 Alltech Associates, Inc. W. R. Grace & Co. is not affiliated with, associated with, endorsed by, or sponsored by Teledyne Isco, Inc. and Biotage AB. All rights reserved. The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn. and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification.

09/2010 M342a

W. R. Grace & Co.-Conn. 7500 Grace Drive, Columbia, MD 21044

Grace Davison Discovery Sciences Regional Headquarters:

IN THE AMERICAS:
2051 Waukegan Road
Deerfield, IL 60015-1899
Tel: 1 847.948.8600
Fax: 1 847.948.1078
Email: discoverysciences@grace.com

IN EUROPE:
Brandstraat 12
B-9160 Lokeren, Belgium
Tel: +32 (0)9-340-65 65
Fax: +32 (0)9-340-65 60
Email: discoverysciences.EU@grace.com

IN ASIA:
19th Floor, K.Wah Center
1010 Huai Hai Zhong Road
Shanghai, 200031 PRC
Tel: 86 21 5467 4678
Fax: 86 21 5405 1500
Email: dsbiz.asia@grace.com

IN INDIA:
17 Commerce Center
Opposite Krishna Hospital, Paud Road
Pune, 411038 India
Tel: +91 20 6644 9900
Fax: +91 20 2544 1740
Email: pune@grace.com

IN AUSTRALIA/NEW ZEALAND:
2 Kerr Court
Rowville, 3178
Victoria, Australia
Tel: +61 3 9237 6100
Fax: +61 3 9237 6101
Email: discoverysciences.AU@grace.com

GRACE