

Column Heaters

Alltech® and Jetstream Models

Jetstream Plus 2



- Peltier column thermostat heats to 85°C and cools to 5°C
- Temperature may be fixed or changed via step or linear gradient
- Large column compartment holds five columns up to 1" o.d. and 35cm in length

Specifications

Access to Columns	Door
Temperature Range:	5–85°C
Column Capacity (Excluding Guards):	5
Maximum Column o.d.:	1"
Maximum Column Length:	35cm
Cavity Dimensions:	15.4" H x 3.3" W x 2" L (39.2cm H x 8.5cm W x 5cm L)
Temperature Accuracy:	±0.2°C
Temperature Stability:	±0.1°C
Operating Voltage:	100–245VAC, 50/60 Hz, 100W
Warranty:	1 year

Ordering Information

Model	Jetstream Plus 2	
	Part No.	
Block heater with one block set and RS232	—	
Single output controller, long column oven and RS232	—	
Dual output controller, 2 long column ovens, and RS232 interface	—	
Heater/Chiller with RS232	—	
Jetstream 2 Plus Heater/Chiller, 75cm cable	531102, dual voltage	
Jetstream 2 Plus Heater/Chiller, 150cm cable	531202, dual voltage	